


COVID-19 Content Pack for Dragon Medical

To assist customers addressing the COVID-19 pandemic, Nuance has built a **content pack** specific to COVID-19 to improve provider efficiency when using Dragon Medical to capture crucial data points for COVID-19 patients. This content pack can be distributed to Dragon Medical One and Dragon Medical Network Edition users through the Nuance Management Center (NMC).

Deployment instructions

1. Download the [Nuance COVID-19 Content Pack XML file](#).
2. Choose a deployment strategy. Nuance recommends deployment across sites for broad sharing. Customers can choose to deploy by site, group, or user.
3. Import the Nuance COVID-19 Content Pack XML file, as shown below.
4. Select the **Speech Recognition** tab.
5. Select the **Manage** icon.
6. Ensure **Auto-texts** is selected.
7. Click the **Import** button.



8. Select **Nuance Management Server format (XML)**, browse to the XML file location and click **Next**.

Import Auto-texts

File to Import

Select an XML file containing the auto-text information:

Dragon format (XML)

Nuance Management Server format (XML)

Nuance COVID-19 Content Pack.xml

9. Click **Select All** to import all Auto-texts and click **Next**.

Import Auto-texts

Auto-texts to Import

The listed auto-texts will be imported.

<input type="checkbox"/>	Auto-text Name	Owner Type	Owner Name	Description
<input checked="" type="checkbox"/>	COVID Actual Exposure			Seconda
<input checked="" type="checkbox"/>	COVID Acute Bronchitis			Seconda
<input checked="" type="checkbox"/>	COVID Acute Respiratory...			Seconda
<input checked="" type="checkbox"/>	COVID ARDS			Seconda
<input checked="" type="checkbox"/>	COVID Assessment			Index Au
<input checked="" type="checkbox"/>	COVID Bronchitis NOS			Seconda
<input checked="" type="checkbox"/>	COVID Lower Respiratory...			Seconda
<input checked="" type="checkbox"/>	COVID Other Respiratory I...			Seconda
<input checked="" type="checkbox"/>	COVID Pneumonia			Seconda

10. If prompted (not all customers will see this), select your organization.


Import Auto-texts

Organization

Select an organization. The auto-texts will be associated with this organization.

Your Organization Name

11. Select your chosen **deployment level**.



12. Click **Next** to import the Auto-texts.

13. Once the import is complete, click **Next** to verify no errors occurred during import.

14. Click **Finish**.

Once imported, the COVID-19 Auto-texts are available for use. In addition to deploying the COVID-19 Content Pack, we recommend customers add “COVID-19” to the Dragon Medical vocabulary using the appropriate instructions below.

- [Instructions for adding COVID-19 to the Dragon Medical vocabulary through Nuance Management Center](#)
- [Instructions for adding COVID-19 to an individual profile vocabulary in Dragon Medical One](#)
- [Instructions for adding COVID-19 to an individual profile vocabulary in Dragon Medical Network Edition](#)

How to use the Content Pack:

- [How to use Auto-texts instructional video](#)
- [Nuance COVID-19 Content Pack for Dragon Medical instructional video](#) provides an overview of content and instructions for use.